Go to the main content Go to menu Go to search



Information and resources for people with disabilities

Disaboom Network: Community | Dating | Jobs

Search search disaboom.com Search

Skip Menu

**Sports Mobility Tech Travel Parenting Education Entertainment Veterans** 

All **Conditions**  Wellness Women's Health

Job Advocacy **Disability Rights** 

**Products & Services** 

**Video Resources** 

Home > Mobility > Adaptive Driving Equipment Information > Electric Wheelchair Carrier Essentials

# **Electric Wheelchair Carrier Essentials**

by Misty Faucheux

text size



email share print



Electric wheelchairs or power chairs tend to be heavier than regular wheelchairs. So, an electric wheelchair carrier for transporting these chairs tends to be a bit more complicated than a wheelchair carrier for a manual chair. For example, a top-mounted wheelchair carrier

wouldn't work for these machines. The weight of an electric wheelchair would probably dent the top of most vehicles.

#### **Electric Wheelchair Carrier Options**

There are, however, carriers made specifically for electric wheelchairs. Before wheelchair owners go out and buy an electric wheelchair carrier, however, they should take their vehicle to a lift expert to find out which type of carrier will work best with their vehicle. Plus, electric wheelchair lifts often require a lot more work than a traditional wheelchair carrier. Since they are motorized, the vehicle's storage space may need to be drilled into in order to wire the device to the car battery.

The most common electric wheelchair carrier for scooters and motorized electric wheelchairs is the type that loads your electric wheelchair into the cargo space or third-seat area of a vehicle. Basically, the unit is attached to the inside of the vehicle and usually latches to the floor.

With the push of a button, the electric wheelchair carrier lifts the power chair up and out the rear of the vehicle, and it will then gently lower the wheelchair to the ground. And, while the car is in motion, it will hold the scooter in place.

#### **Electric Wheelchair Ramp**

If wheelchair owners have room, they can also consider using a <u>wheelchair ramp</u>. For most scooters, a 10 foot ramp would be sufficient in order to drive the wheelchair into a vehicle. But, asking the advice of a professional would probably be best before purchasing the ramp.

Once the wheelchair is in or out of the vehicle, the ramp will fold up into the rear of the vehicle. And, then the power wheelchair can be secured within the car or truck via a variety of methods. Ramps work just as well for electric wheelchair lifts as those that mount to the car itself. Further, ramps can be removed when they are not in use to provide additional space within the vehicle.

Photo Courtsey of **E-Z Carrier** 

Back to Adaptive Driving Equipment Information

### email print



Receive our newsletter

Newsletter Signup enter email address Sign Up

Related Content



Top 8 Used Wheelchair Vans Sites



# Defensive Driving with Adaptive Equipment



### Wheelchair Carriers: What You Need to Know



Back on the Road Again: Adaptive Vehicles for SCIs



Advertisement

WellPoint.

Neglistered Trackmark, Diversity Inc Media LLC.

Registered Trackmark, Diversity Inc Media LLC.

Popular topics

### **Adaptive Sports**

Adaptive sports is an exceptionally popular topic among Disaboom readers. Whether looking for information about how to participate, where to... <u>View More</u>

## Cancer

Cancer comes in many forms, and has just as many possible approaches to treatment. If you or someone you love has received a cancer... <u>View More</u>

### **Emergency Preparedness**

Emergency preparedness for people with disabilities is critically important, as has been demonstrated in recent natural disasters. Do you... <u>View More</u>

Our Company

About Us

Contact Us

Advertise With Us

Terms of Use

Our Network

**DisaboomJobs** 

**DisaboomLive** 

Lovebyrd

Resources

**Disability Organizations** 

**Disability Scholarships** 

Outreach

Assistive Technology Devices

What is Disaboom?

Disaboom is the leading resource for disability information and real-life articles about people with disabilities. Our broad range of topics, including health conditions, lifestyle, and helpful resources, help you create the life you want.